粉铸 (Powercast)

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# PowerCast



A carbon-free, fine-grain phosphate bonded investment for casting crown and bridge and ceramic alloys. This versatile all-purpose investment has been specially formulated to use with a rapid burnout technique. It can be used both with a metal ring and a ringless technique.

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1. Pattern Preparation:

- Attach sprue to the wax pattern using pliable sticky wax.
- Mount patterns on crucible former base.
- Spray with Whip Mix SMOOTHEx Wax Pattern Cleaner. Gently blow excess cleaner from pattern.

2. Preparation:

- Prepare liquid at recommended concentration as required following the chart below.
- If metal ring is being used, line the ring with a ceramic or paper liner positioned 6 mm short on each end.
- Tack liner in place with pliable sticky wax.
- Soak lined ring in water for a minimum of one minute, then shake well (5 times). Lined ring can be soaked in SMOOTHEx for more expansion.
- Place ring onto former base.

3. Mixing:

Suggested Liquid/Powder Ratio:

60 g	90 g	100 g
14 mL	21 mL	23 mL

- Rinse bowl out with water — shake out excess.
- Add measured liquid into bowl then add weighed powder.
- Hand mix well for 20 seconds.
- Mechanical mix under vacuum, slow speed (350–450 RPM) for 90–120 seconds.

4. Benchset:

- If using ringless technique, remove mold from ring and base after 20–30 minutes. Then benchset for an additional 25 minutes.
- NOTE: Minimum benchset time is 45 minutes.
- If using metal ring, remove former base from ring prior to placing mold in oven.
- Trim glaze off top of mold and rinse under tap water.

5. Burnout:

Option “A” — Rapid Technique if patterns and sprues are wax:

- Place molds in preheated oven at 700–850° C (1300–1550° F). (If higher temperature is warranted, place mold in preheated oven at 370 °C (700 °F) and then raise to final temperature.)
- Heat soak at final temperature for 30 minutes, add 10 minutes per additional mold.
- For rapid burnout of single ringless molds, place in preheated oven at 430°C (800°F) and then raise to final temperature.

Option “B” — Two Stage Technique if patterns or sprues contain plastic:

- Place mold in furnace at room temperature and raise temperature at a rate of 8° C (15° F) per minute to 430° C (800° F).
- Hold for 30 minutes.
- Raise to maximum temperature at a rate of 14° C (25° F) per minute.
- Hold for 30 minutes.

6. Casting:

- Cast alloy with oxygen/propane torch or induction machine according to alloy manufacturer’s instructions.
- Suggested turns on centrifugal casting machine:  
Crown and bridge gold alloy.....1-2 turns  
Ceramic gold alloy .....2-3 turns  
High palladium and non-precious .....3 turns

Note:

- Recommended Liquid: Special Liquid Concentrate Plus
- Use separate mixing bowls for phosphate and gypsum investments.
- For best results, keep the liquid and powder at room temperature of about 15–24 °C (60–75 °F).
- When using bulk, reseal container to avoid fast set and breakdown.
- Special Liquid Concentrate Plus may freeze during shipment in cold weather. If liquid should have crystallized upon arrival, it will no longer be usable. Keep liquid from freezing.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer).

Physical Properties:

Liquid*/Powder Ratio	23 mL /100 g
Working Time	8–9 minutes
Setting Expansion	0.8%
Thermal Expansion	1.0%
Compressive Strength, Wet	700 psi (5 MPa)
*Suggested concentration of Special Liquid Concentrate Plus is 80% (4 parts liquid to 1 part water).	

Special Liquid Concentrate Plus Recommended Concentrations:

Alloy	Ringless	Metal Ring w/ceramic liner*	Metal Ring w/paper liner	Burnout Temperature
Chrome-Cobalt	80–100%	100%	100%	900°C (1650°F)
Nickel-Beryllium	75–100%	90–100%	100%	845°C (1550°F)
Palladium (2–8% Gold)	50–60%	50–60%	60–75%	815°C (1500°F)
Silver-Palladium	60–75%	75–80%	80–100%	790°C (1450°F)
Ceramic Gold (≥45%)	50–75%	50–75%	75–80%	760°C (1400°F)
Crown & Bridge Gold	60–75%	60–75%	75–80%	650°C (1200°F)

\* Recommended concentrations are approximate and can be adjusted to optimize fit. To increase expansion, use more liquid and less water. To decrease expansion, use less liquid and more water. Always maintain total liquid / water volume.

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